



♥ PEDIATRIC

Assessing or diagnosing a metabolic disorder commonly requires several tests. Global Metabolomic Assisted Pathway Screen, commonly known as Global MAPS, is a unifying test for analyzing hundreds of metabolites to identify changes or irregularities in biochemical pathways. Let Global MAPS guide you to an answer.

BAYLOR
GENETICS


GLOBAL MAPS™
Global Metabolomic
Assisted Pathway Screen

Diagnose a broad range of
metabolic disorders with a
single test, Global MAPS

Global MAPS is a large scale, semi-quantitative metabolomic profiling screen that analyzes disruptions in both individual analytes and pathways related to biochemical abnormalities.

Using state-of-the-art technologies, Global Metabolomic Assisted Pathway Screen (Global MAPS) provides small molecule metabolic profiling to identify >700 metabolites in human plasma, urine, or cerebrospinal fluid. Global MAPS identifies inborn errors of metabolism (IEMs) that would ordinarily require many different tests. This test defines biochemical pathway errors not currently detected by routine clinical or genetic testing.

IEMs are inherited metabolic disorders that prevent the body from converting one chemical compound to another or from transporting a compound in or out of a cell.



These processes are necessary for essentially all bodily functions. Most IEMs are caused by defects in the enzymes that help process nutrients, which result in an accumulation of toxic substances or a deficiency of substances needed for normal body function. Making a swift, accurate diagnosis of an IEM and prescribing the appropriate diet or medication is critical in preventing brain damage, organ damage, and even death.





Global MAPS Indications For Testing

AUTISM SPECTRUM DISORDER

DEVELOPMENTAL DELAY

VARIANTS OF UNCERTAIN CLINICAL
SIGNIFICANCE IN A GENE KNOWN TO BE
INVOLVED IN SMALL MOLECULE METABOLISM

FAILURE TO THRIVE

HYPOGLYCEMIA

HYPOTONIA

NON-SYNDROMIC INTELLECTUAL DISABILITY

RECURRENT VOMITING

SEIZURES

SPEECH/LANGUAGE DELAY

UNDIFFERENTIATED PHENOTYPE
POSSIBLY RELATED TO PERTURBATION
IN A BIOCHEMICAL PATHWAY



Global MAPS provides patients assist in the diagnosis of a recognize new IEMs never

The possibility of discovery by using Global MAPS is considerable, including the identification of new inborn errors of metabolism and regulatory factors for metabolic genes, as well as previously unknown metabolic associations/ disruptions with known disorders. Global MAPS offers a broad range of analyses in a single metabolic screen, requiring less sample volume from the patient making it more convenient and cost-effective.

Global MAPS is a unique broad screening test that can detect disorders involving metabolism of amino acids, organic acids, fatty acid oxidation, vitamin cofactors, pyrimidine biosynthesis, creatine biosynthesis, and urea cycle metabolism, among other known disorders.

Metabolites range in size and include, but are not limited to:

AMINO ACIDS


FATTY ACIDS

ORGANIC ACIDS

NEUROTRANSMITTERS

NUCLEOTIDES




BILE ACIDS



ents with a single test that can broad range of disorders and before described.

Testing Options

GLOBAL METABOLOMIC ASSISTED PATHWAY SCREEN (GLOBAL MAPS)

Test Code	4900	4901	4902
Specimen Type			
TAT (Days)	21	21	21

Specimen Requirements & Shipping Conditions



Send 1-2 cc of plasma. Draw blood in an EDTA (purple top) tube(s) and separate plasma as soon as possible, freezing immediately. Store the specimen frozen at -20°C. Specimen may be stored frozen up to 7 days.

Ship frozen sample in insulated container, with 3 -5 lbs. dry ice, by overnight courier.



Send 1-2 cc of cerebrospinal fluid. Store the specimen frozen at -20°C. Specimen may be stored frozen for up to 7 days.

Ship frozen sample in insulated container, with 3 -5 lbs. dry ice, by overnight courier.



Send 3-5 cc of a random urine. Do not add preservatives. Store the specimen frozen at -20°C.

Ship frozen sample in insulated container, with 3 -5 lbs. dry ice, by overnight courier.



Plasma



Urine



Cerebrospinal Fluid

EDTA (Purple Top)




40 YEARS OF INNOVATION



4 MILLION+ CLINICAL TESTS PERFORMED



1 MILLION+ FAMILIES HELPED



3 THOUSAND+ TESTS OFFERED



1 MISSION IMPROVE HEALTHCARE THROUGH GENETICS

**Baylor Genetics pioneered the history of genetic testing.
Now, we're leading the way in precision medicine.**

Baylor Genetics is a joint venture of Miraca Holdings, Inc. and Baylor College of Medicine, including the #1NIH funded Department of Molecular and Human Genetics. A pioneer of precision medicine for nearly 40 years, Baylor Genetics now offers a full spectrum of clinically relevant genetic testing, access to world-renowned experts, and the confidence to provide patients with the best care.

1.800.411.4363
BAYLORGENETICS.COM



BAYLOR
GENETICS

Global MAPS Disorders List





PLASMA (CONFIRMED DISORDERS DETECTED BY GLOBAL MAPS)

Urea Cycle Disorders

Argininemia
Argininosuccinic aciduria
Citrullinemia
Ornithine transcarbamylase deficiency
Orotic aciduria

Organic Acidemias

2-hydroxyglutaric acidemia (likely L-form)
3-hydroxyisobutyryl-CoA hydrolase deficiency (HIBCH)
3-hydroxy-3-methylglutaryl(HMG)-CoA lyase deficiency
3-methylcrotonyl-CoA carboxylase (3-MCC) deficiency
Cobalamin biosynthesis disorders
Ethylmalonic encephalopathy
Glutaric acidemia type I
Holocarboxylase synthetase deficiency
Isovaleric acidemia
Methylmalonic acidemia
Propionic Acidemia

Fatty Acid Oxidation Disorders

Short chain acyl-CoA dehydrogenase (SCAD) deficiency
Medium chain acyl-CoA dehydrogenase (MCAD) deficiency
Very long chain acyl-CoA dehydrogenase (VLCAD) deficiency

Amino Acid Disorders

Citrin deficiency
Classical homocystinuria (cystathionine β-synthetase deficiency)
Glycine encephalopathies
Hyperphenylalaninemia
Lysinuric protein intolerance
Maple syrup urine disease
Phenylketonuria
Serine biosynthesis disorders (phosphoserine aminotransferase deficiency, phosphoglycerate dehydrogenase deficiency)
Tyrosinemia type I

Other

Adenylosuccinate lyase deficiency
Aromatic L-amino acid decarboxylase deficiency
β-Ureidopropionase deficiency
Citrate transporter deficiency
Creatine biosynthesis defects (GAMT and AGAT deficiencies)
GABA transaminase deficiency
Galactosemia
Glycerol kinase deficiency
Hyperornithinemia-hyperammonemia-homocitrullinuria (HHH) syndrome
MTHFR deficiency
Mitochondrial neurogastrointestinal encephalopathy (MNGIE)
Peroxisome biogenesis disorders/Zellweger spectrum disorders
Primary carnitine deficiency
Pyridoxine-dependent epilepsy
Ribose-5-phosphate isomerase deficiency
Smith-Lemli-Opitz syndrome
Spondyloepimetaphyseal dysplasia, Genevieve type
Thiamine transporter deficiency
Trimethyllysine hydroxylase epsilon deficiency
Transaldolase deficiency
Transketolase deficiency
Urocanase deficiency (benign condition)



PLASMA (DISORDERS THAT SHOULD BE DETECTED BUT HAVE NOT BEEN ASSESSED)

For this group, we routinely detect one or more plasma analytes that are well established biomarkers for disease.

2-Methylbutyryl-CoA Dehydrogenase Deficiency
3-Hydroxyacyl-CoA Dehydrogenase (SCHAD) Deficiency
AMACR Deficiency
Beta-Ketothiolase Deficiency
Canavan Disease
Carbamoyl Phosphate Synthetase I Deficiency
Carnitine-Acylcarnitine Translocase Deficiency
Carnitine Palmitoyltransferase I (CPT1) Deficiency
Carnitine Palmitoyltransferase II (CPT2) Deficiency
Combined Malonic And Methylmalonic Aciduria

Dihydropyrimidinase Deficiency
Dimethylglycine Dehydrogenase Deficiency
Fructose-1,6-Bisphosphatase Deficiency
Glutaric Acidemia II
Gyrate Atrophy Of Choroid And Retina
Hereditary Fructose Intolerance
Holocarboxylase Synthetase Deficiency
Hypermethioninemia due to S-Adenosylhomocysteine Hydrolase Deficiency
Hypermethioninemia Due To Adenosine Kinase Deficiency
Hyperoxaluria Type I
Hyperoxaluria Type II

Hyperprolinemia, Type I
Hyperprolinemia, Type II
Lathosterolosis
Lesch-Nyhan Syndrome
Malonyl-CoA Decarboxylase Deficiency
Molybdenum Cofactor Deficiency
N-Acetylglutamate Synthase (NAGS) Deficiency
Phosphoribosylpyrophosphate Synthetase (PRPPS) Superactivity
Phosphoserine Phosphatase Deficiency
Purine Nucleoside Phosphorylase Deficiency
Succinyl-CoA:3-Oxoacid CoA Transferase (SCOT) Deficiency



PLASMA (DISORDERS THAT MAY BE DETECTED BUT HAVE NOT BEEN VALIDATED)

We identify one or more compounds that are predicted to be relevant to the disorder on the basis of their position within the affected metabolic pathway.

α-Aminoacidic Aciduria	Congenital Bile Acid Synthesis Defect 1 (CBAS1)	Hyperphenylalaninemia, BH4-Deficient, D (PCBD Deficiency)
2-Methyl-3-Hydroxybutyryl-CoA Dehydrogenase Deficiency	Congenital Bile Acid Synthesis Defect 3 (CBAS3)	Hyperuricemic Nephropathy, Familial Juvenile
3-Methylglutaconic Aciduria, Type I	Dihydropyrimidine Dehydrogenase Deficiency	Intrinsic Factor Deficiency; IFD
3-Methylglutaconic Aciduria, Type III	Essential Fructosuria (considered benign)	Kelley-Seegmiller Syndrome
3-Methylglutaconic Aciduria, Type V	Familial Hypercholanemia (Bile Acid Biosynthesis Disorder)	Methionine Adenosyltransferase Deficiency
3-Methylglutaconic Aciduria Type VI (with deafness, encephalopathy, and Leigh-like syndrome)	Glutamine Deficiency, Congenital	Prolidase Deficiency
D-2-Hydroxyglutaric Aciduria Type 1	Glycine N-Methyltransferase Deficiency	Succinic Semialdehyde Dehydrogenase Deficiency
D-2-Hydroxyglutaric Aciduria Type 2	Hydroxykynureninuria	Transcobalamin II Deficiency
Adenine Phosphoribosyltransferase (APRT) Deficiency	Hyperphenylalaninemia, BH4-Deficient, A (PTS Deficiency)	Tyrosinemia, Type II
Asparagine Synthetase Deficiency	Hyperphenylalaninemia, BH4-Deficient, B (GTP Cyclohydrolase Deficiency)	Tyrosinemia, Type III
Carnosinemia	Hyperphenylalaninemia, BH4-Deficient, C (DHPR Deficiency)	Xanthinuria, Type I

Limitations to testing: Disorders we presume we cannot detect on this small molecule test

- Congenital disorders of protein glycosylation
- Glycogen storage diseases
- Oligosaccharidoses
- Mucopolysaccharidoses
- Other lysosomal storage diseases



URINE (THE FOLLOWING DISORDERS DETECTABLE ONLY BY URINE)

- Creatine transporter deficiency (detectable in urine Global MAPS only)



CEREBROSPINAL FLUID (THE FOLLOWING DISORDERS DETECTABLE ONLY BY CEREBROSPINAL FLUID)

- NFU1 deficiency (detectable in CSF Global MAPS only)