Before You Start

Do not eat, drink (including water), smoke, or chew gum at least 30 minutes before collecting your sample.

Collecting a sample takes about 2-5 minutes.

Sample must be shipped via FedEx® within 24 hours of collection.

Scan the QR code to watch a quick video that shows each step for collecting and sending your sample.

Instructions are available in other languages. Please see the manufacturer instructions included in your kit for more details.

Kit Contents




(2) ORAcollect®•Dx Collection Kits each Containing a Buccal Swab Tube with Sponge Tip

(1) Manufacturer Instructions


ORAcollect·Dx

© DNA Genotek Inc.
All rights reserved

(1) Español
Instrucciones de Usuario



(4) Tube Labels



(1) Biohazard Bag with Absorbent Pad



(1) Prepaid FedEx® Return Mailer



If the kit is missing any items or if they appear damaged, please contact Client Services at 1.800.411.4363.

BAYLOR GENETICS

How to Collect & Ship a Buccal Sample for Testing



WHAT YOU RECEIVED

You have received an individual kit or a family collection kit to collect samples. Each person being tested should have their own kit.


If you are collecting samples from more than one family member, put each person's sample in its own box so they do not get mixed up. The family collection kit packaging sleeve may be thrown away.



If you have any questions or need assistance, please contact Client Services at 1.800.411.4363.


Sample Collection

1



Open package and remove collector without touching sponge tip. Place sponge as far back in the mouth as comfortable and rub along the lower gums (see close up image) in a back and forth motion. Gently rub the gums 10 times. If possible, avoid rubbing the teeth.

2



Gently repeat rubbing motion on the opposite side of the mouth along the lower gums for an additional 10 times.

3




Hold the tube upright to prevent the liquid inside the tube from spilling. Unscrew the blue cap from the collection tube without touching the sponge.

Illustrations and instructions appearing in 'Sample Collection' (steps 1–5) are used with permission from © DNA Genotek Inc. All rights reserved.


Sample Collection

4



Turn the cap upside down, insert the sponge into the tube and close cap tightly.

5




Invert the capped tube and shake vigorously 10 times.

Repeat steps 1-5 for each buccal swab tube.

Labeling

6



Using the blank labels provided and clearly print:

- Name of Person's Sample Being Collected
- Date of Birth
- Date of Collection

Date of Collection is the date the sample was obtained. Use one label per specimen container and attach it to the container(s).

If information is missing, or the tube label does not match what your provider sent, we cannot run the test.


Packaging & Shipping

7



Place the specimen container(s) into the biohazard bag that contains an absorbent pad. Make sure the biohazard bag is completely sealed. Each sample goes in its own biohazard bag. If included, fold the test order form and insert it into the slip pouch, located on the exterior, clear side of the bag.

8



Place the sealed biohazard bag into the collection kit box. Place the closed box in the FedEx® mailer and seal it. Multiple boxes can be put into the same mailer. Unused return mailers can be thrown away.

9

FedEx®

SHIP WITHIN 24 HRS OF COLLECTION (MON-FRI)

Drop off at a FedEx® location or call 800-GoFedEx to schedule a pickup. No account is required; use the tracking number from the FedEx label.